

In Vivo CNS Receptor Occupancy & In vitro Receptor Autoradiography



- **Beta-Imager**


- Can detect all beta emitting isotopes (^3H , ^{14}C , ^{32}P , ^{33}P , ^{35}S , ^{125}I)
- Shortens the time needed for assessment of tritiated ligands imaging time from weeks to hours when comparing to film
- Therefore is deal for:
 - **Receptor occupancy studies:** an assay that could provide crucially important preclinical data for preclinical and clinical candidate selection, as a screening tool.
 - Ex Vivo: Only cold compound is administered (PO, IP, SC)
 - In vivo: Both radiolabeld (IV) and cold compound (PO, IP, SC) is administered
 - **Receptor autoradiography:** Determine the tissue distributions of target of interest. To obtain ED50, Ki, Kd, and Bmax values.



Redefining Drug Discovery Through Innovation

Quantitative Autoradiography

There Is A 33% Decrease in D2 Receptor Expression in 9 Month Old Heterozygous Q175KI Male Mice Model of HD




Dopamine D2 receptor saturation binding in the striatum of 9 month old z_Q175KI HET mice and their WT littermates. Calculated Bmax values for WT and z_Q175KI HET animals shows a loss in D2 receptor expression at 9 months of age (33%, **** $p < 0.0001$).




Redefining Drug Discovery Through Innovation

Ex Vivo Occupancy

Measurement of PDE10A inhibitor Target Occupancy



Vehicle



1mg/kg

3mg/kg

10mg/kg


Ex vivo binding: PDE10 inhibitor, P.O., 1hr




Ex Vivo RO/Efficacy Relationship

Sazetidine ($\alpha 4\beta 2$ nAChR partial agonist) After ip Dosing in Mouse


Behavior: Forced Swim



Ex Vivo Binding



Saz (mg/kg, i.p, 30min)



High correspondence between behaviorally active doses of sazetidine in the forced swim test and $\alpha 4\beta 2$ nicotinic acetylcholine receptor (nAChR) occupancy in the brain. Specifically, sazetidine was inactive in the forced swim test at 0.3 mg/kg (i.p) and showed very low receptor occupancy at this dose (<10%). In contrast, doses of sazetidine (1 and 3 mg/kg, i.p.) that produced robust behavioral activity in the forced swim test, corresponded with high levels of occupancy at the $\alpha 4\beta 2$ nAChRs (~40% and ~90% respectively).



Caldarone et al, Psychopharmacology (Berl). 2011



Redefining Drug Discovery Through Innovation

In vivo Occupancy

Brain Receptor Occupancy of CTP-354, a Selective Benzodiazepine Receptor ($GABA_A$) Modulator in Rats



Tracer: 3H -Ro15-1788 (Flumazenil) administered by tail vein injection.




Redefining Drug Discovery Through Innovation

^{35}S GTP γ S Binding
Assessment of Functional Activation of G_i/G_o -Coupled Receptors

μ -Opioid Receptors agonist DAMGO stimulates $^{35}\text{SGTP}\gamma\text{S}$ binding in the rat spinal cord, amygdala, and striatum/accumbens.

- The $^{35}\text{SGTP}\gamma\text{S}$ assay allows for the assessment of functional activation of Gi/Go-coupled receptors by measuring bound $^{35}\text{SGTP}\gamma\text{S}$ in brain sections using autoradiography.
- This assay is useful for investigating biased agonism in m-opioid receptors (MORs).







Redefining Drug Discovery Through Innovation

^{14}C -2-Deoxyglucose (^{14}C -2-DG) Autoradiography
Glucose Utilization as a Function of Neurogenic Activity

Measuring Glucose Consumption After Acute i.p. Administration of 20 mg/kg Clozapine Using ^{14}C -2-DG Autoradiography

- The ^{14}C -2-DG autoradiographic method can be used for measuring glucose utilization as a function of neurogenic activity.
- This assay is useful for identifying brain regions activated by specific drug agents, and therefore, can act as a surrogate for functional MRI.



Fifteen minutes after clozapine dosing, rats were i.v. tail vein injected with 25 $\mu\text{Ci}/\text{kg}$ of ^{14}C -2-DG. Rats were sacrificed 45 minutes post i.v. injection of ^{14}C -2-DG.